



IAL-SYSTEM **DUO**

CROSS-LINKED HYALURONIC ACID + ACP,™ 2.5%

NEW



The 4th Dimension in Aesthetic Medicine



BIOREVOLUMETRY® 



PHITOGEN HOLDING

VOLUME AND TISSUE RE-STRUCTURING

BIOREVOLUMETRY[®] green injections

IAL-SYSTEM[®] DUO

CROSS-LINKED HYALURONIC ACID + ACP[™], 2.5%

IAL-SYSTEM DUO is a bio-revolumetric product that belongs to the category of injection treatments aimed at complete rejuvenation of the face through volumetric correction and tissue re-structuring.


It is an innovative, apyrogenic, clear and transparent, highly viscous gel, which is the result of a physical association of two gels obtained through the cross-linking of Hyaluronic Acid produced by fermentation. Its prolonged permanence in site offers a long-lasting filling effect with a slow release of free Hyaluronic Acid, which – thanks to its natural hydrating and stimulating properties, improves skin turgor and elasticity and maintains them over a long period of time.

Restores volume and re-structures tissues

LONG-LASTING BIO-REVOLUMETRIC ACTION IN LOCALIZED HYPO-VOLUMES

IAL-SYSTEM DUO is indicated in case of localized hypo-volumes (LH) which can be mainly found in more or less young patients, and are basically characterized by a poor presence of the deep and/or superficial adipose layer.

CONCENTRATION	25mg/ml
TECHNOLOGY	Mu.C.H. [™] (Multi-cross-linked Hyaluronan)
COMPOSITION	75% cross-linked Hyaluronic Acid (with BDDE) 25% ACP (auto-cross-linked Hyaluronic Acid without BDDE)
DYNAMIC VISCOSITY	on average 100 Pas
OSMOLARITY	330 mOsm/L
RESIDUAL BDDE	< 2 ppm
ENDOTOXINS	< 2,4 EU/ml
STORAGE	2-8°C



CE 0373 IAL-SYSTEM[®] DUO
CROSS-LINKED HYALURONIC ACID + ACP[™], 2.5%



Manufacturer
FIDIA FARMACEUTICI Spa
www.fidiapharma.com



PHITOGEN HOLDING Spa
www.phitogenspa.eu